

LUCAS™ Chest Compression System

MYTH BUSTED! LUCAS *does* fit large patients!



Weight: 180 lb
Height: 6 ft



Weight: 330 lb
Height: 6 ft, 5 in



Weight: 320 lb
Height: 5 ft, 10 in





Proper Placement is on Upper Chest

Do not assume the LUCAS Chest Compression System can't fit larger patients. Many large patients carry their weight through their mid-section. When LUCAS is properly applied, it rests over the upper chest or sternum, with enough room for the compressor to give 1.5 to 2 inch compressions. As long as the unit snaps into the backboard and the suction cup does not compress the patient's chest in the start position, it will operate as intended.

LUCAS 1

CHEST COMPRESSION SYSTEM



Patients Eligible for Treatment:

- Sternum height of 7.5 – 11.9 inches (19 – 30.3 cm)
- Maximum sternum width of 17.7 inches (45 cm)

LUCAS 2

CHEST COMPRESSION SYSTEM



Patients Eligible for Treatment:

- Sternum height of 6.7–11.9 inches (17 – 30.3 cm)
- Maximum chest width: 17.7 inches (45 cm)

The use of LUCAS is not restricted by patient weight.

LUCAS is designed and manufactured by Jolife AB in Lund, Sweden and distributed worldwide by Physio-Control, Inc.

For further information, please contact your sales representative at 800.442.1142 or visit our website at www.physio-control.com



Physio-Control, Inc.
11811 Willows Road NE
P. O. Box 97006
Redmond, WA 98073-9706 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com